

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Pongo abelii Lysosome-associated membrane glycoprotein 5 (LAMP5), partial BYT-ORB1096472

Artikelname	Recombinant Pongo abelii Lysosome-associated membrane glycoprotein 5 (LAMP5), partial
Artikelnummer	BYT-ORB1096472
Hersteller Artikelnummer	orb1096472
Alternativnummer	BYT-ORB1096472-1,BYT-ORB1096472-100,BYT-ORB1096472-20
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	This Recombinant Pongo abelii Lysosome-associated membrane glycoprotein 5 (LAMP5), partial spans the amino acid sequence from region 30-234aa. Purity: Greater than 90% as determined by SDS-PAGE....
Molekulargewicht	30.4 kDa
UniProt	Q5R5V2
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Reinheit	Greater than 90% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	EQEVENLSGLSTNPEKDIFVVRENGTTCLMAEFAAKFIVPYDVWASNYVDLITE QADIALTRGAEVGRCGHSESELQVFWVDLAYALKMLFVKESHNMSKGPEET WRLSKVQFVYDSSEKTHFKDAVSAGKHTANSHHLSALVTPAGKSYECQAQQ TISLASSDPQKVTMILSAVHIQPFDIISDFVFSEEHKCAVDEREQLEE
Anwendungsbeschreibung	Biological Origin: Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference