

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant *Saccharomyces cerevisiae* Dolichyl-phosphate-mannose--protein mannosyltransferase 2 (PMT2), partial BYT-ORB1096509

Artikelname	Recombinant <i>Saccharomyces cerevisiae</i> Dolichyl-phosphate-mannose--protein mannosyltransferase 2 (PMT2), partial
Artikelnummer	BYT-ORB1096509
Hersteller Artikelnummer	orb1096509
Alternativnummer	BYT-ORB1096509-1,BYT-ORB1096509-100,BYT-ORB1096509-20
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	This Recombinant <i>Saccharomyces cerevisiae</i> Dolichyl-phosphate-mannose--protein mannosyltransferase 2 (PMT2), partial spans the amino acid sequence from region 311-609aa. Purity: Greater than 85% as determined by SDS-PAGE....
Molekulargewicht	38.0 kDa
UniProt	<a href="#">P31382</a>
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Bakers yeast)
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	LLSHSGTGDANMPSLFQARLVGSDVGQGPRDIALGSSVVSIGNQALGGSLH SHIQTYPDGSNQQQVTCYGYKDANNEWFFNRERGLPSWSENETDIEYLKPGT SYRLVHKSTGRNLHHPVAAPVSKTQWEVSGYGDNVVGDNKDNWVIEIMDQ RGDEDPEKLHTLTTSFRIKNLEMGCYLAQTGNSLPEWGFRRQEVVCMKNPFK RDKRTWWNIETHENERLPPRPEDFQYPKTNFLKDFIHLNLMAMMATNNA
Anwendungsbeschreibung	Biological Origin: <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Bakers yeast). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference