

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human coronavirus OC43 Non-structural protein 2a (2a) BYT-ORB1096667

|                          |  |
|--------------------------|--|
| Artikelname              | Recombinant Human coronavirus OC43 Non-structural protein 2a (2a)  |
| Artikelnummer            | BYT-ORB1096667   |
| Hersteller Artikelnummer | orb1096667   |
| Alternativnummer         | BYT-ORB1096667-1,BYT-ORB1096667-100,BYT-ORB1096667-20  |
| Hersteller               | Biorbyt  |
| Kategorie                | Proteine/Peptide   |
| Produktbeschreibung      | This Recombinant Human coronavirus OC43 Non-structural protein 2a (2a) spans the amino acid sequence from region 1-278aa. Purity: Greater than 90% as determined by SDS-PAGE....   |
| Molekulargewicht         | 33.3 kDa   |
| UniProt                  | <a href="#">Q80872</a>   |
| Puffer                   | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Quelle                   | Human coronavirus OC43 (HCoV-OC43)   |
| Reinheit                 | Greater than 90% as determined by SDS-PAGE.  |
| Formulierung             | Liquid or Lyophilized powder   |

|                        |  |
|------------------------|--|
| Sequenz                | MAVAYADKPNHFINFPLTHFQGFVLNYKGLQFQILDEGVDCIKIQTAPHISLTML<br>DIQPEDYKSVDVAIQEVIDDMHWGDGFGQIKFENPHILGRCIVLDVKGVEELHD<br>DLVNYIRDKGCVADQSRKWIGHCTIAQLTDAALSIKENVDFINSMQFNKITIN<br>PSSPARLEIVKLGAEEKDGFYETIVSHWMGIRFEYTSPTDKLAMIMGYCCLDVV<br>RKELEEGDLPENDDDAWFKLSYHYENNSWFFRHVYRKSFH   |
| Anwendungsbeschreibung | Biological Origin: Human coronavirus OC43 (HCoV-OC43). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference |