

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Sheep anti Rabbit Myokinase, Clone: [Polyclonal], Monoclonal NMB-NE106/7S

| | |
|--------------------------|--|
| Artikelname | Sheep anti Rabbit Myokinase, Clone: [Polyclonal], Monoclonal |
| Artikelnummer | NMB-NE106/7S |
| Hersteller Artikelnummer | NE106/7S |
| Alternativnummer | NMB-NE106/7S |
| Hersteller | NordicMubio |
| Wirt | Sheep |
| Kategorie | Antikörper |
| Produktbeschreibung | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and... |
| Klonalität | Monoclonal |
| Klon-Bezeichnung | [Polyclonal] |
| UniProt | P00569 |
| Puffer | Purified hyperimmune sheep IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile distilled |
| Quelle | Myokinase isolated and purified from rabbit muscle. Freund's complete adjuvant is used in the first step of the immunization procedure. |

| | |
|------------------------|---|
| Formel | Purified hyperimmune sheep IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Anwendungsbeschreibung | Indirect immunofluorescence,ELISA,Dot blot,Immunoblotting. |