

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISpot Plus: Human IFN-gamma (ALP) MAB-3420-4AST-10

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ELISpot Plus: Human IFN-gamma (ALP)                                                                                                                                                                                                                           |
| Artikelnummer            | MAB-3420-4AST-10                                                                                                                                                                                                                                              |
| Hersteller Artikelnummer | 3420-4AST-10                                                                                                                                                                                                                                                  |
| Alternativnummer         | MAB-3420-4AST-10                                                                                                                                                                                                                                              |
| Hersteller               | MABTECH                                                                                                                                                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISpot                                                                                                                                                                                                                                                       |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Konjugation              | ALP                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | ELISpot Plus for enumeration of cells secreting IFN-gamma, This kit is ideal for users who want a convenient and sensitive assay. The assay is designed for the enumeration of cells secreting human IFN-gamma. The kit includes 96-well ELISpot strip pla... |
| Detektionsbereich        | Streptavidin-ALP                                                                                                                                                                                                                                              |
| NCBI                     | <a href="#">3458</a>                                                                                                                                                                                                                                          |
| Pubmed                   | <a href="https://pubmed.ncbi.nlm.nih.gov/3131436">https://pubmed.ncbi.nlm.nih.gov/3131436</a> target=_blank rel=nofollow>3131436                                                                                                                              |
| Target-Kategorie         | IFN-gamma                                                                                                                                                                                                                                                     |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>Complementary Products:</p> <p><a href="https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0">https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefta-cd4-human-3617-1-0">https://www.mabtech.com/products/peppool-cefta-cd4-human-3617-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-human-3627-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-human-3627-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0</a>. Contents: 7-B6-1, biotin,Streptavidin-ALP enzyme conjugate for ELISpot,BCIP/NBT-plus substrate for ELISpot,Anti-CD3 mAb (CD3-2),Pre-coated ELISpot plate, strip (mAb 1-D1K)</p> |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|