

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISA Flex: Human IL-2 (HRP) MAB-3445-1H-6

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ELISA Flex: Human IL-2 (HRP)                                                                                                                                                                                                                                  |
| Artikelnummer            | MAB-3445-1H-6                                                                                                                                                                                                                                                 |
| Hersteller Artikelnummer | 3445-1H-6                                                                                                                                                                                                                                                     |
| Alternativnummer         | MAB-3445-1H-6                                                                                                                                                                                                                                                 |
| Hersteller               | MABTECH                                                                                                                                                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISA                                                                                                                                                                                                                                                         |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | ELISA Flex kit for quantification of IL-2 in solution, The human IL-2 ELISA Flex kit is ideal for users who want the flexibility of setting up their own ELISA assay. This product enables quantification of human IL-2 in solution, e.g. in cell culture ... |
| Detektionsbereich        | 4-400 pg/ml                                                                                                                                                                                                                                                   |
| NCBI                     | <a href="https://pubmed.ncbi.nlm.nih.gov/29291917">3558</a>                                                                                                                                                                                                   |
| Pubmed                   | <a href="https://pubmed.ncbi.nlm.nih.gov/29291917">29291917</a>                                                                                                                                                                                               |
| Target-Kategorie         | IL-2                                                                                                                                                                                                                                                          |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | Complementary Products: <a href="https://www.mabtech.com/products/elisa-diluent-3652-d2-0">https://www.mabtech.com/products/elisa-diluent-3652-d2-0</a> , <a href="https://www.mabtech.com/products/elisa-substrate-tmb-hrp-3652-f10-0">https://www.mabtech.com/products/elisa-substrate-tmb-hrp-3652-f10-0</a> . Contents: MT2A91/2C95,MT8G10, biotin ,Streptavidin-ALP enzyme conjugate for ELISpot,Recombinant human IL-2 ELISA standard,Standard reconstitution buffer |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|