

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISpot Flex: Mouse IL-17A (ALP) MAB-3521-2A

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ELISpot Flex: Mouse IL-17A (ALP)                                                                                                                                                                                                                              |
| Artikelnummer            | MAB-3521-2A                                                                                                                                                                                                                                                   |
| Hersteller Artikelnummer | 3521-2A                                                                                                                                                                                                                                                       |
| Alternativnummer         | MAB-3521-2A                                                                                                                                                                                                                                                   |
| Hersteller               | MABTECH                                                                                                                                                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISpot                                                                                                                                                                                                                                                       |
| Spezies Reaktivität      | Mouse                                                                                                                                                                                                                                                         |
| Konjugation              | ALP                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | Basic ELISpot kit for enumeration of cells secreting IL-17A, The mouse IL-17A ELISpot kit is ideal for users who want the flexibility of setting up their own assay. The kit enables enumeration of cells secreting mouse IL-17A, and contains a matched p... |
| Detektionsbereich        | Streptavidin-ALP                                                                                                                                                                                                                                              |
| NCBI                     | <a href="#">16171</a>                                                                                                                                                                                                                                         |
| Pubmed                   | <a href="https://pubmed.ncbi.nlm.nih.gov/26263201" rel="nofollow" target="_blank">26263201</a>                                                                                                                                                                |
| Target-Kategorie         | IL-17A                                                                                                                                                                                                                                                        |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | Complementary Products:<br><a href="https://www.mabtech.com/products/elispot-substrate-bcipnbt-plus-alp-3650-10-0">https://www.mabtech.com/products/elispot-substrate-bcipnbt-plus-alp-3650-10-0</a> , <a href="https://www.mabtech.com/products/elispot-plates-white-3654-wp">https://www.mabtech.com/products/elispot-plates-white-3654-wp</a> . Contents: IL17-I,MT2270, biotin,Streptavidin-ALP enzyme conjugate for ELISpot |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|