

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISpot Flex: Mouse IgA (ALP) MAB-3865-2A

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ELISpot Flex: Mouse IgA (ALP)                                                                                                                                                                                                                                 |
| Artikelnummer            | MAB-3865-2A                                                                                                                                                                                                                                                   |
| Hersteller Artikelnummer | 3865-2A                                                                                                                                                                                                                                                       |
| Alternativnummer         | MAB-3865-2A                                                                                                                                                                                                                                                   |
| Hersteller               | MABTECH                                                                                                                                                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISpot                                                                                                                                                                                                                                                       |
| Spezies Reaktivität      | Mouse                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | Basic ELISpot for enumeration of IgA secreting B cells, Mouse IgA ELISpot Flex enables enumeration of B cells secreting total mouse IgA as well as antigen-specific mouse IgA. It contains matched quantities of monoclonal capture and detection antibody... |
| Detektionsbereich        | Streptavidin-ALP                                                                                                                                                                                                                                              |
| NCBI                     | <a href="#">238447</a>                                                                                                                                                                                                                                        |
| Pubmed                   | <a href="https://pubmed.ncbi.nlm.nih.gov/24958702">https://pubmed.ncbi.nlm.nih.gov/24958702</a> target=_blank rel=nofollow>24958702                                                                                                                           |
| Target-Kategorie         | IgA                                                                                                                                                                                                                                                           |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | Complementary Products:<br><a href="https://www.mabtech.com/products/elispot-plates-white-3654-wp">https://www.mabtech.com/products/elispot-plates-white-3654-wp</a> ,<br><a href="https://www.mabtech.com/products/elispot-substrate-bciphbt-plus-alp-3650-10-0">https://www.mabtech.com/products/elispot-substrate-bciphbt-plus-alp-3650-10-0</a> . Contents: MT45A,MT39A, biotin,Streptavidin-ALP<br>enzyme conjugate for ELISpot,R848,Recombinant mouse IL-2 |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|