

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

SBP | Sedoheptulose-1,7-bis phosphatase positive control/quantitation standard **AGR-AS15-2873S**

Artikelname	SBP Sedoheptulose-1,7-bis phosphatase positive control/quantitation standard
Artikelnummer	AGR-AS15-2873S
Hersteller Artikelnummer	AS15-2873S
Alternativnummer	AGR-AS15-2873S
Hersteller	Agrisera
Kategorie	Sonstiges
Applikation	WB
Produktbeschreibung	SBPase (Sedoheptulose-1,7-bis phosphatase) is a chloroplast enzyme involved in the carbon reduction of the Calvin cycle, part of carbohydrate metabolism.Source of SBPase standard: SBPase standard source is Sphingomonas wittichii strain RW1, overexpre...
Molekulargewicht	27 kDa
Formulierung	Lyophilized
Antibody Type	Secondary Antibody
Application Verdünnung	Standard curve: 3 loads are recommended (0.5, 2 and 4 µl).For most applications a sample load of 0.2 µg of chlorophyll/well will give a RbcL signal in this range.Positive control: a 2 µl load per well is optimal for most chemiluminescent detection systems

Anwendungsbeschreibung	Concentration: after re-constitution with sterile milliQ water final concentration of the standard is 0.15 pmoles/μl Protein standard buffer composition: Glycerol 10%, Tris Base 141 mM, Tris HCl 106 mM, LDS 2%, EDTA 0.51 mM, SERVA Blue G250 0.22 mM, Phenol Red 0.175 mM, pH 8.5, 0.1mg/ml PefaBloc protease inhibitor (Roche), 50 mM DTT. This standard is ready-to-load and does not require any additions or heating. It needs to be fully thawed and thoroughly mixed prior to using. Avoid vigorous vortexing, as buffers contain detergent. Following mixing, briefly pulse in a microcentrifuge to collect material from cap. This standard is stabilized and ready and does not require heating before loading on the gel. Please note that this product contains 10% glycerol and might appear as liquid but is provided lyophilized. Allow the product several minutes to solubilize after adding water. Mix thoroughly but gently Take extra care to mix thoroughly before each use, as the proteins tend to settle with the more dense layer after freezing.
------------------------	--