

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit IgG anti-Chicken IgG (H+L)-Biotin, MinX none DNA-SEC-182802

Artikelname	Rabbit IgG anti-Chicken IgG (H+L)-Biotin, MinX none
Artikelnummer	DNA-SEC-182802
Hersteller Artikelnummer	SEC-182802
Alternativnummer	DNA-SEC-182802
Hersteller	dianova
Wirt	Rabbit
Kategorie	Antikörper
Applikation	ELISA,IHC,WB
Spezies Reaktivität	Gallus
Immunogen	Anti-Chicken IgY / IgG whole molecule was produced by repeated immunization with Chicken IgY / IgG whole molecule in rabbit.
Konjugation	Biotin
Format	IgG
Spezifität	IgG (H+L)
Minimale Kreuzreaktivität (MinX)	no cross-adsorbtion
Produktbeschreibung	Anti-Chicken IgG Biotin Antibody generated in rabbit detects chicken IgY. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as w...

Klonalität	Polyclonal
Konzentration	1.5 mg/mL
Isotyp	Ig
Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Biotin Conjugated Anti-Chicken IgG / IgY secondary antibody was prepared from monospecific antiserum by immunoaffinity chromatography using Chicken IgG / IgY coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin, anti-Rabbit Serum and Chicken IgG / IgY. Some light chain cross reactivity may be observed against chicken immunoglobulins.
Formulierung	Lyophilized
Formel	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% NaN3
Target-Kategorie	Chicken
Antibody Type	Secondary Antibody
Application Verdünnung	ELISA Dilution: 1:500,000, Immunohistochemistry Dilution: 1:1,000 - 1:5,000, Western Blot Dilution: 1:2,000 - 1:20,000
Anwendungsbeschreibung	Secondary antibodies are ideal for western blotting, Immunohistochemistry, ELISA, Fluorescence Microscopy, Flow Cytometry as well as other antibody detection methods.