

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### 3T3 cells, Unconjugated NMB-X1005

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | 3T3 cells, Unconjugated                                                                                                                                                                                                                                       |
| Biozol Catalog Number      | NMB-X1005                                                                                                                                                                                                                                                     |
| Supplier Catalog Number    | X1005                                                                                                                                                                                                                                                         |
| Alternative Catalog Number | NMB-X1005                                                                                                                                                                                                                                                     |
| Manufacturer               | NordicMubio                                                                                                                                                                                                                                                   |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                              |
| Application                | WB                                                                                                                                                                                                                                                            |
| Conjugation                | Unconjugated                                                                                                                                                                                                                                                  |
| Product Description        | Total cell lysates are useful as both positive and negative controls in immunoblotting. A431, a human epidermoid carcinoma cell line, and the EGF-stimulated A431 lysates are used as negative and positive controls, respectively, when studying the phos... |
| Concentration              | See vial for concentration                                                                                                                                                                                                                                    |
| Buffer                     | Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.                                                                                                                             |
| Form                       | Unconjugated                                                                                                                                                                                                                                                  |
| Formula                    | Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.                                                                                                                             |
| Application Notes          | Suggested use: A volume of 10-25ul is suggested for immunoblotting.                                                                                                                                                                                           |