

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

A431 cells stimulated with EGF, Unconjugated NMB-X1003

Article Name	A431 cells stimulated with EGF, Unconjugated
Biozol Catalog Number	NMB-X1003
Supplier Catalog Number	X1003
Alternative Catalog Number	NMB-X1003
Manufacturer	NordicMubio
Category	Proteine/Peptide
Application	WB
Conjugation	Unconjugated
Product Description	A431 cells are a positive control for EGFR as they express high levels of EGFR. They also used for cell cycle studies and can be stimulated with EGF to stimulate EGFR and are therefore a very useful model to study EGFR downstream signalling cascades ...
Concentration	See vial for concentration
Buffer	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Form	Unconjugated
Formula	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Application Notes	Suggested use: A volume of 10-25ul is suggested for immunoblotting.