

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Chicken anti Sheep IgG (H+L), conjugated with FITC, Fluorescein, Gallus NMB-ASIG-F

|                            |                                                                         |
|----------------------------|-------------------------------------------------------------------------|
| Article Name               | Chicken anti Sheep IgG (H+L), conjugated with FITC, Fluorescein, Gallus |
| Biozol Catalog Number      | NMB-ASIG-F                                                              |
| Supplier Catalog Number    | ASIG-F                                                                  |
| Alternative Catalog Number | NMB-ASIG-F                                                              |
| Manufacturer               | NordicMubio                                                             |
| Host                       | Gallus                                                                  |
| Category                   | Antikörper                                                              |
| Application                | IHC                                                                     |
| Species Reactivity         | Sheep                                                                   |
| Conjugation                | Fluorescein                                                             |
| Format                     | IgY                                                                     |
| Target Specificity         | IgG (H+L)                                                               |
| Cross-Adsorption (MinX)    | no cross-adsorbtion                                                     |
| Product Description        | Chicken anti Sheep IgG (H+L), conjugated with FITC, Gallus 0,5 mg...    |
| Concentration              | See vial for concentration.                                             |
| Buffer                     | Phosphate Buffered Saline, 50% Glycerol with 0.075% sodium azide        |

|                   |                                                                                                                                                                                                                                                                  |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Source            | Chickens were immunized with highly purified sheep IgG isolated from pooled normal sheep serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid p |
| Purity            | Egg-yolk derived IgY Affinity Purified                                                                                                                                                                                                                           |
| Form              | Fluorescein conjugated                                                                                                                                                                                                                                           |
| Formula           | Phosphate Buffered Saline, 50% Glycerol with 0.075% sodium azide                                                                                                                                                                                                 |
| Antibody Type     | Secondary Antibody                                                                                                                                                                                                                                               |
| Application Notes | Optimal working dilutions should be determined for your particular assay conditions.                                                                                                                                                                             |