

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Rabbit IgG Fab, Clone: [Polyclonal], Monoclonal NMB-GAR/FAB

Article Name	Goat anti Rabbit IgG Fab, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-GAR/FAB
Supplier Catalog Number	GAR/Fab
Alternative Catalog Number	NMB-GAR/FAB
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Species Reactivity	Rabbit
Conjugation	Unconjugated
Format	Antiserum
Target Specificity	IgG (F(ab') ₂)
Cross-Adsorption (MinX)	no cross-adsorbtion
Product Description	The antiserum is reacting with the Fab subunit of intact IgG, IgA and IgM and other Ig classes of both light chain types, with their Fab or F(ab) ₂ subunits and with free light chains of kappa and lambda type as tested in immunoelectrophoresis and dou...
Clonality	Monoclonal
Clone Designation	[Polyclonal]

Buffer	Delipidated, heat inactivated lyophilized stable whole serum No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added. It is reconstituted by adding 1 ml st
Source	Purified Fab of normal rabbit IgG. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formula	Delipidated, heat inactivated lyophilized stable whole serum No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added.
Antibody Type	Secondary Antibody
Application Notes	Precipitation assays. In immunoelectrophoresis use 2 μ l or equivalent against 120 μ l antiserum. In double radial immunodiffusion use a rosette arrangement with 10 μ l antiserum in a 3 mm diameter centre well and 2 μ l serum samples (neat and diluted) in 2 mm diameter peripheral wells.