

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse peroxidase-anti-peroxidase, Clone: [Polyclonal], Monoclonal NMB-M/PAP

| | |
|----------------------------|---|
| Article Name | Mouse peroxidase-anti-peroxidase, Clone: [Polyclonal], Monoclonal |
| Biozol Catalog Number | NMB-M/PAP |
| Supplier Catalog Number | M/PAP |
| Alternative Catalog Number | NMB-M/PAP |
| Manufacturer | NordicMubio |
| Host | Mouse |
| Category | Antikörper |
| Product Description | Peroxidase staining of unlabelled mouse IgG antibodies and tissue sections.... |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| Buffer | Hyperimmune mouse IgG against horseradish peroxidase reacted with peroxidase lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divid |
| Source | Highly purified horseradish peroxidase. |
| Formula | Hyperimmune mouse IgG against horseradish peroxidase reacted with peroxidase lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). |