

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Rabbit anti Enterobacter aerogenes Glycerol dehydrogenase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE068/BIO**

Article Name	Rabbit anti Enterobacter aerogenes Glycerol dehydrogenase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE068/BIO
Supplier Catalog Number	NE068/Bio
Alternative Catalog Number	NMB-NE068/BIO
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	Biotin
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...
Clonality	Monoclonal
Clone Designation	[Polyclonal]
UniProt	<a href="#">G0EAP7</a>
Buffer	Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile distilled

Source	Glycerol dehydrogenase isolated and purified from <i>Enterobacter aerogenes</i> . Freund's complete adjuvant is used in the first step of the immunization procedure.
Form	Biotin
Formula	Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.