

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rabbit anti Porcine Acylase, Clone: [Polyclonal], Monoclonal NMB-NE163/PAB

|                            |   |
|----------------------------|---|
| Article Name               | Rabbit anti Porcine Acylase, Clone: [Polyclonal], Monoclonal  |
| Biozol Catalog Number      | NMB-NE163/PAB   |
| Supplier Catalog Number    | NE163/Pab   |
| Alternative Catalog Number | NMB-NE163/PAB   |
| Manufacturer               | NordicMubio   |
| Host                       | Rabbit  |
| Category                   | Antikörper  |
| Product Description        | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...   |
| Clonality                  | Monoclonal  |
| Clone Designation          | [Polyclonal]  |
| Buffer                     | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitute |
| Source                     | Acylase I isolated and purified from porcine kidney. Freund's complete adjuvant is used in the first step of the immunization procedure.  |

|                   |  |
|-------------------|--|
| Formula           | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Application Notes | Immunoprecipitation.   |