

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rabbit anti Porcine Glutamic-pyruvic transaminase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE064/BIO

Article Name	Rabbit anti Porcine Glutamic-pyruvic transaminase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE064/BIO
Supplier Catalog Number	NE064/Bio
Alternative Catalog Number	NMB-NE064/BIO
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	Biotin
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...
Clonality	Monoclonal
Clone Designation	[Polyclonal]
UniProt	<a href="#">P13191</a>
Buffer	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile dis

Source	Glutamic-pyruvic transaminase isolated and purified from porcine heart. Freund's complete adjuvant is used in the first step of the immunization procedure.
Form	Biotin
Formula	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.