

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Gamma-Aminobutyric Acid Receptor Subunit Alpha-2 (GABRA2) ELISA KBS-KBH0487

|                            |                                                                                                                                                                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Gamma-Aminobutyric Acid Receptor Subunit Alpha-2 (GABRA2) ELISA                                                                                                                                                                                |
| Biozol Catalog Number      | KBS-KBH0487                                                                                                                                                                                                                                                  |
| Supplier Catalog Number    | KBH0487                                                                                                                                                                                                                                                      |
| Alternative Catalog Number | KBS-KBH0487-ELISA,KBS-KBH0487-1                                                                                                                                                                                                                              |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                          |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                  |
| Species Reactivity         | Human                                                                                                                                                                                                                                                        |
| Product Description        | The GENLISA Human Gamma-Aminobutyric Acid Receptor Subunit Alpha-2 (GABRA2) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Preci... |
| Range                      | 0.156 - 10 ng/ml                                                                                                                                                                                                                                             |
| Sensitivity                | 0.05 ng/ml                                                                                                                                                                                                                                                   |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                           |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                   |