

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### 2 \* Rapid Taq Master Mix VZB-P222-02

|                            |                                                                                                                                                                    |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | 2 * Rapid Taq Master Mix                                                                                                                                           |
| Biozol Catalog Number      | VZB-P222-02                                                                                                                                                        |
| Supplier Catalog Number    | P222-02                                                                                                                                                            |
| Alternative Catalog Number | VZB-P222-02-15                                                                                                                                                     |
| Manufacturer               | Vazyme Biotech                                                                                                                                                     |
| Category                   | Sonstiges                                                                                                                                                          |
| Product Description        | 1. The PCR extension speed can be set at 15 sec/kb and even 1 sec/kb for DNA segments within 1kb. 2. The master mix remains stable after 50 freeze-thaw cycles.... |
| Application Notes          | Rapid PCR                                                                                                                                                          |