

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

LPS from Salmonella typhimurium S-form (TLRGRADE) (Ready-to-Use) ENZ-ALX-581-011-L002

Article Name	LPS from Salmonella typhimurium S-form (TLRGRADE) (Ready-to-Use)
Biozol Catalog Number	ENZ-ALX-581-011-L002
Supplier Catalog Number	ALX-581-011-L002
Alternative Catalog Number	ENZ-ALX-581-011-L002-2
Manufacturer	Enzo Life Sciences
Category	Sonstiges
Source	Smooth (S)-form LPS, isolated and purified from Salmonella typhimurium by modification of the phenol water extraction and PCP method, converted to the uniform salt form and dissolved in sterile pyrogen-free double distilled water.
Purity	Absence of detectable protein or DNA contaminants with agonistic TLR activity.
Form	Liquid. Sterile, ready-to-use solution in pyrogen-free double distilled water.
Target	TLR4