

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

LPS from Salmonella abortus equi S-form (TLRGRADE) (Ready-to-Use) (Biotin) **ENZ-ALX-581-150-R500**

Article Name	LPS from Salmonella abortus equi S-form (TLRGRADE) (Ready-to-Use) (Biotin)
Biozol Catalog Number	ENZ-ALX-581-150-R500
Supplier Catalog Number	ALX-581-150-R500
Alternative Catalog Number	ENZ-ALX-581-150-R500-500
Manufacturer	Enzo Life Sciences
Category	Sonstiges
Conjugation	Biotin
Product Description	Useful reagent for evaluating LPS binding and cell activation in white blood cells....
Source	Smooth (S)-form LPS, isolated and purified from Salmonella abortus equi by modification of the phenol water extraction and PCP method. LPS-Biotin was prepared using the biotin reagent biotinamidocaproate N-hydroxysuccinimide ester. Briefly, LPS at 10mg/m
Purity	Absence of detectable protein or DNA contaminants with agonistic TLR activity.
Form	Liquid. Sterile, ready-to-use solution in pyrogen-free double-distilled water.
Target	TLR4