

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Matrix metalloproteinase (MMP) multipack-2 ENZ-BML-AK014-0001

Article Name	Matrix metalloproteinase (MMP) multipack-2
Biozol Catalog Number	ENZ-BML-AK014-0001
Supplier Catalog Number	BML-AK014-0001
Alternative Catalog Number	ENZ-BML-AK014-0001-1
Manufacturer	Enzo Life Sciences
Category	Proteine/Peptide
Application	FA, HTS
Product Description	Matrix metalloproteinase (MMP) multipack-2...
Form	Liquid. Each enzyme is supplied in 50mM TRIS, pH 7.5, containing 5mM CaCl ₂ , 300mM NaCl, 20μM ZnCl ₂ , 30% glycerol, 0.5% Brij-35.
Target	MMP MMP MMP-3 MMP-7 MMP-10 MMP-11 MMP-12
Application Notes	Enzyme regulation and kinetics, comparative studies of substrate or inhibitor specificities, cleavage of target proteins. Assay or digest conditions can vary widely, but concentrations for the MMP enzymes can range between 10 and 300nM, or higher. Reacti