

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SUMO-2 (human), (recombinant) (His-tag) **ENZ-BML-UW9205-0500**

Article Name	SUMO-2 (human), (recombinant) (His-tag)
Biozol Catalog Number	ENZ-BML-UW9205-0500
Supplier Catalog Number	BML-UW9205-0500
Alternative Catalog Number	ENZ-BML-UW9205-0500-500
Manufacturer	Enzo Life Sciences
Category	Proteine/Peptide
Product Description	SUMO-2 (human), (recombinant) (His-tag)...
Molecular Weight	15.4 kDa
UniProt	P61956
Source	Produced in E. coli. Human SUMO-2 is fused to an N-terminal His-tag.
Form	Liquid. In 20mM HEPES, pH 8.0, containing 50mM NaCl, 1mM DTT.
Target	SUMO
Application Notes	Uses: SUMO-2 modification of specific proteins in vitro.Demonstrate novel proteins are potential targets for SUMOylation under in vitro conditions.Generate substrates for deSUMOylating enzymes.