

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Fluorescein Recombinant monoclonal antibody (4-4-20 enhanced) (Mouse Fab fragment kappa) (His-tag) ENZ-ENZ-ABS493-0200**

|                            |                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------|
| Article Name               | Fluorescein Recombinant monoclonal antibody (4-4-20 enhanced) (Mouse Fab fragment kappa) (His-tag)    |
| Biozol Catalog Number      | ENZ-ENZ-ABS493-0200                                                                                   |
| Supplier Catalog Number    | ENZ-ABS493-0200                                                                                       |
| Alternative Catalog Number | ENZ-ENZ-ABS493-0200-200                                                                               |
| Manufacturer               | Enzo Life Sciences                                                                                    |
| Host                       | Mouse                                                                                                 |
| Category                   | Antikörper                                                                                            |
| Application                | FC, WB                                                                                                |
| Immunogen                  | Fluorescein I-bovine serum albumin (FII-BSA).                                                         |
| Product Description        | Fluorescein Recombinant monoclonal antibody (4-4-20 enhanced) (Mouse Fab fragment kappa) (His-tag)... |
| Clonality                  | Monoclonal                                                                                            |
| Clone Designation          | [4-4-20 enhanced]                                                                                     |
| Purity                     | Immobilized metal affinity chromatography purified.                                                   |
| Form                       | Liquid. In PBS containing 0.02% Proclin 300.                                                          |
| Target                     | Fluorescein                                                                                           |
| Application Notes          | Negative control                                                                                      |