

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Fluorescent NADP/NADPH Detection Kit ICT-NADPH100-3

|                            |                                                                                                                                                                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Fluorescent NADP/NADPH Detection Kit                                                                                                                                                                                                                |
| Biozol Catalog Number      | ICT-NADPH100-3                                                                                                                                                                                                                                      |
| Supplier Catalog Number    | NADPH100-3                                                                                                                                                                                                                                          |
| Alternative Catalog Number | ICT-NADPH100-3                                                                                                                                                                                                                                      |
| Manufacturer               | ImmunoChemistry Technologies                                                                                                                                                                                                                        |
| Category                   | Kits/Assays                                                                                                                                                                                                                                         |
| Product Description        | A non-fluorescent detection reagent is reduced in the presence of NADPH to produce its fluorescent analog and NADP. NADP is further converted to NADPH via an enzyme-coupled reaction. The kit specifically detects NADP/NADPH and not NAD/NADH.... |