

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Napsin A (Lung Adenocarcinoma Marker) Antibody, Clone: [NAPSA/1865R], Rabbit, Monoclonal NBT-9476-RBM5-P0

|                            |  |
|----------------------------|--|
| Article Name               | Recombinant Napsin A (Lung Adenocarcinoma Marker) Antibody, Clone: [NAPSA/1865R], Rabbit, Monoclonal   |
| Biozol Catalog Number      | NBT-9476-RBM5-P0   |
| Supplier Catalog Number    | 9476-RBM5-P0   |
| Alternative Catalog Number | NBT-9476-RBM5-P0-20,NBT-9476-RBM5-P0-100   |
| Manufacturer               | NeoBiotechnologies   |
| Host                       | Rabbit   |
| Category                   | Antikörper   |
| Application                | IHC, WB  |
| Species Reactivity         | Human, Rabbit, Rat   |
| Immunogen                  | Recombinant human Napsin-A protein fragment (aa189-299) (exact sequence is proprietary)  |
| Product Description        | Napsin is a pepsin-like aspartic proteinase connected with maturation of surfactant protein B. There are two closely related napsins, napsin A and napsin B. Napsin A is expressed as a single chain protein. Immunohistochemical studies revealed high exp... |
| Clonality                  | Monoclonal   |
| Clone Designation          | [NAPSA/1865R]  |
| Molecular Weight           | 37kDa  |
| NCBI                       | <a href="#">9476</a>   |

|                   |   |
|-------------------|---|
| UniProt           | <a href="#">O96009</a>  |
| Form              | 200ug/ml of recombinant MAb Purified by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.   |
| Antibody Type     | Recombinant Monoclonal Antibody   |
| Application Notes | Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT)(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes), Optimal dilution f |