

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Collagen I for plate coating - Bovine Skin ROC-001-G70-MF6

|                            |                                                                                                                                    |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Collagen I for plate coating - Bovine Skin                                                                                         |
| Biozol Catalog Number      | ROC-001-G70-MF6                                                                                                                    |
| Supplier Catalog Number    | 001-G70-MF6                                                                                                                        |
| Alternative Catalog Number | ROC-001-G70-MF6                                                                                                                    |
| Manufacturer               | Rockland Immunochemicals                                                                                                           |
| Category                   | Kits/Assays                                                                                                                        |
| Application                | SDS-PAGE                                                                                                                           |
| Conjugation                | Unconjugated                                                                                                                       |
| Concentration              | 1.0 mg/mL by dry weight                                                                                                            |
| Buffer                     | 0.013M HCl                                                                                                                         |
| Form                       | Liquid (sterile filtered)                                                                                                          |
| Application Notes          | Bovine Collagen Type I solution has been tested in SDS-page and is suitable for plate coating. Enough to coat 10 (ten) T75 plates. |