

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Bovine Serum Albumin 30% Solution ROC-BSA-30-0050

Article Name	Bovine Serum Albumin 30% Solution
Biozol Catalog Number	ROC-BSA-30-0050
Supplier Catalog Number	BSA-30-0050
Alternative Catalog Number	ROC-BSA-30-0050
Manufacturer	Rockland Immunochemicals
Category	Zellen/Zellkultur
Conjugation	Unconjugated
Concentration	30 % by dry weight
Buffer	0.85% (w/v) Sodium Chloride
Form	Liquid (sterile filtered)
Application Dilute	ELISA: User Optimized, FLISA: User Optimized, IHC: User Optimized, IF Microscopy: User Optimized, IP: User Optimized, WB: User Optimized
Application Notes	BOVINE SERUM ALBUMIN 30% Solution is suitable for use in protease sensitive assays such as RIA, EIA and nucleic acid hybridization, use as a stabilizing agent for proteins and enzymes, including dilute solutions of antibody, and use as a blocking agent t