

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Chemiluminescent FemtoMax(TM) Super Sensitive HRP Substrate ROC-FEMTOMAX-110

Article Name	Chemiluminescent FemtoMax(TM) Super Sensitive HRP Substrate
Biozol Catalog Number	ROC-FEMTOMAX-110
Supplier Catalog Number	FEMTOMAX-110
Alternative Catalog Number	ROC-FEMTOMAX-110
Manufacturer	Rockland Immunochemicals
Category	Sonstiges
Application	WB
Conjugation	Unconjugated
Concentration	1X
Form	Liquid - clear, colorless, odorless
Application Dilute	ELISA: 1X, WB: 1X
Application Notes	Prepare FemtoMax(TM) Super Sensitive Chemiluminescent HRP Substrate for use in microwell or membrane applications by mixing 1 mL of Luminol Reagent (Reagent A) with 1 mL of Reaction Buffer (Reagent B). Mix well. Protect from light. Larger or smaller volumes