

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### SARS-CoV-2 (COVID-19) Helicase Recombinant Protein PRS-92-765

|                            |                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------|
| Article Name               | SARS-CoV-2 (COVID-19) Helicase Recombinant Protein                                         |
| Biozol Catalog Number      | PRS-92-765                                                                                 |
| Supplier Catalog Number    | 92-765                                                                                     |
| Alternative Catalog Number | PRS-92-765-0.05                                                                            |
| Manufacturer               | ProSci                                                                                     |
| Category                   | Proteine/Peptide                                                                           |
| NCBI                       | <a href="#">43740578</a>                                                                   |
| UniProt                    | <a href="#">P0DTC1</a>                                                                     |
| Buffer                     | Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.5. |
| Source                     | E. coli                                                                                    |
| Purity                     | Greater than 80% as determined by reducing SDS-PAGE.                                       |
| Form                       | Liquid                                                                                     |
| Application Notes          | This recombinant protein can be used for biological assays. For research use only.         |