

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Thymus Membrane Tumor Lysate PRS-XBL-10984

|                            |                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------|
| Article Name               | Thymus Membrane Tumor Lysate                                                                                        |
| Biozol Catalog Number      | PRS-XBL-10984                                                                                                       |
| Supplier Catalog Number    | XBL-10984                                                                                                           |
| Alternative Catalog Number | PRS-XBL-10984-0.1                                                                                                   |
| Manufacturer               | ProSci                                                                                                              |
| Category                   | Zellen/Zellkultur                                                                                                   |
| Application                | WB                                                                                                                  |
| Concentration              | 2 mg/mL or better                                                                                                   |
| Buffer                     | HEPES pH 7.9, MgCl <sub>2</sub> , KCl, EDTA, Sucrose, Glycerol, Sodium deoxycholate, and NP-40.                     |
| Form                       | Liquid                                                                                                              |
| Application Notes          | Thymus membrane protein lysate is for use in Western blotting, 10 µg to 20 µg per lane is recommended for mini gel. |