

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant *Oryza sativa* subsp. *japonica* Mitogen-activated protein kinase 5 (MPK5), Unconjugated, *E. coli* BIM-RPC22956

Article Name	Recombinant <i>Oryza sativa</i> subsp. <i>japonica</i> Mitogen-activated protein kinase 5 (MPK5), Unconjugated, <i>E. coli</i>
Biozol Catalog Number	BIM-RPC22956
Supplier Catalog Number	RPC22956
Alternative Catalog Number	BIM-RPC22956-20UG,BIM-RPC22956-100UG,BIM-RPC22956-1MG
Manufacturer	Biomatik Corporation
Host	<i>E. coli</i>
Category	Proteine/Peptide
Conjugation	Unconjugated
Product Description	Recombinant <i>Oryza sativa</i> subsp. <i>japonica</i> Mitogen-activated protein kinase 5 (MPK5) is a purified Recombinant Protein. Purity: >85% as determined by SDS-PAGE. Host: <i>E. coli</i> . Endotoxin Level: Not Tested. Species: <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice). Ta...
Molecular Weight	48kDa
Tag	N-Terminal 10Xhis-Tagged And C-Terminal Myc-Tagged
Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Purity	>85% by SDS-PAGE

Sequence	MDGAPVAEFRPTMTHGGRYLLYDIFGNKFEVTNKYQPPIMPIGRGAYGIVCSV MNFETREMVAIKKIANAFNNDMDAKRTLREIKLLRHLDHENIIGIRDVIPPIPQ AFNDVYIATELMDTDLHHIIRSNQELSEHCQYFLYQILRGLKYIHSANVIHRDL KPSNLLLNNANCDLKICDFGLARPSSESMMTEYVVTRWYRAPELLLNSTDYSA AIDVWSVGCIFMELINRQPLFPGRDHMHQMRLITEVIGT
Target	MPK5