

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

RNA Marker Low Easy ABN-R0002

Article Name	RNA Marker Low Easy
Biozol Catalog Number	ABN-R0002
Supplier Catalog Number	R0002
Alternative Catalog Number	ABN-R0002-125
Manufacturer	Abnova
Category	Molekularbiologie
Application	EP
Product Description	The RNA Marker Low Easy is supplied in a ready-to-use mixture of loading dye (containing formamide, EDTA sodium salt, bromphenol blue) and RNAs. It is prepared for denaturing polyacrylamide gel electrophoresis but not agarose gel electrophoresis. The ...
Application Dilute	5 uL is recommended for loading to a well (0.1 ug of each RNA / 5 uL)