

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EpCAM / CD326 (Epithelial Marker) (EGP40/2041R), 1mg/mL, Clone: [EGP40/2041R], Rabbit, Monoclonal BOT-BNUM2041-50

Article Name	EpCAM / CD326 (Epithelial Marker) (EGP40/2041R), 1mg/mL, Clone: [EGP40/2041R], Rabbit, Monoclonal
Biozol Catalog Number	BOT-BNUM2041-50
Supplier Catalog Number	BNUM2041-50
Alternative Catalog Number	BOT-BNUM2041-50-50UL
Manufacturer	Biotium
Host	Rabbit
Category	Antikörper
Application	IHC, WB
Species Reactivity	Human
Immunogen	Recombinant human EpCAM fragment from extracellular domain (around aa100-224) (exact sequence is proprietary)
Product Description	Epitope of this MAb is mapped between aa 202-212 of EGP40, which is a 40-43 kDa transmembrane epithelial glycoprotein, also identified as epithelial specific antigen (ESA), or epithelial cellular adhesion molecule (Ep-CAM). It is expressed on baso-la...
Clonality	Monoclonal
Concentration	1 mg/mL
Clone Designation	[EGP40/2041R]
Molecular Weight	40-43 kDa

UniProt	P16422
Buffer	PBS, no BSA, no azide
Source	Animal
Application Notes	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistology (formalin): 0.5-1 ug/mL for 30 minutes at RT Staining of formalin-fixed tissues requires boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user