

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Canine TM (Thrombomodulin) ELISA Kit EKL-ECA27292

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Canine TM (Thrombomodulin) ELISA Kit                                                                                                                                                                                                                          |
| Biozol Catalog Number      | EKL-ECA27292                                                                                                                                                                                                                                                  |
| Supplier Catalog Number    | ECa27292                                                                                                                                                                                                                                                      |
| Alternative Catalog Number | EKL-ECA27292-96T                                                                                                                                                                                                                                              |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Canine                                                                                                                                                                                                                                                        |
| Product Description        | This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Canine TM. Standards or samples are added to the micro ELISA plate wells and combined with the specific anti... |
| Range                      | 15.62-1000pg/mL                                                                                                                                                                                                                                               |
| Sensitivity                | 5.8pg/mL                                                                                                                                                                                                                                                      |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | TM                                                                                                                                                                                                                                                            |