

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Mouse D2D (D-Dimer) ELISA Kit EKL-EM20506

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Mouse D2D (D-Dimer) ELISA Kit                                                                                                                                                                                                                                 |
| Biozol Catalog Number      | EKL-EM20506                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | EM20506                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | EKL-EM20506-96T                                                                                                                                                                                                                                               |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Mouse                                                                                                                                                                                                                                                         |
| Product Description        | This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Mouse D2D. During the reaction, Mouse D2D in the sample or standard competes with a fixed amount of Mouse D2D on the solid phase ... |
| Range                      | 78.13-5000ng/mL                                                                                                                                                                                                                                               |
| Sensitivity                | 46.88ng/mL                                                                                                                                                                                                                                                    |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | D2D                                                                                                                                                                                                                                                           |