

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **RNase H II, Glycerol-Free SCX-RNH01-E21C-10**

|                            |                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | RNase H II, Glycerol-Free                                                                                                                                           |
| Biozol Catalog Number      | SCX-RNH01-E21C-10                                                                                                                                                   |
| Supplier Catalog Number    | RNH01-E21C-10                                                                                                                                                       |
| Alternative Catalog Number | SCX-RNH01-E21C-10                                                                                                                                                   |
| Manufacturer               | SignalChem Diagnostics                                                                                                                                              |
| Category                   | Sonstiges                                                                                                                                                           |
| Product Description        | RNase H II, Glycerol-Free recognizes DNA/RNA duplexes and cleaves the RNA strand from the 5 end. Its optimal activity is at 70-75C, still active at 50C and 75C.... |