

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### PIP Array - Lipid-Protein Interaction Assay ECH-P-6100

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | PIP Array - Lipid-Protein Interaction Assay                                                                                                                                                                                                                   |
| Biozol Catalog Number      | ECH-P-6100                                                                                                                                                                                                                                                    |
| Supplier Catalog Number    | P-6100                                                                                                                                                                                                                                                        |
| Alternative Catalog Number | ECH-P-6100-2PK,ECH-P-6100-5PK                                                                                                                                                                                                                                 |
| Manufacturer               | Echelon Biosciences                                                                                                                                                                                                                                           |
| Category                   | Molekularbiologie                                                                                                                                                                                                                                             |
| Product Description        | PIP Arrays are 4 x 5 cm hydrophobic membranes that have been spotted with a concentration gradient of all eight phosphoinositides. They are used to determine the relative degree of binding of your protein of interest towards one or more of these lipi... |
| Application Notes          | Categories: Bioassays & Kits. Related: Array, Phosphoinositides. Ship 1 : Shipping Temp. Ship 2: Ambient Temperature                                                                                                                                          |