

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Oculomedin (OCLM) ELISA KBS-KBH10150

|                            |                                                                                                                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Oculomedin (OCLM) ELISA                                                                                                                                                                                                                               |
| Biozol Catalog Number      | KBS-KBH10150                                                                                                                                                                                                                                                        |
| Supplier Catalog Number    | KBH10150                                                                                                                                                                                                                                                            |
| Alternative Catalog Number | KBS-KBH10150-1                                                                                                                                                                                                                                                      |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                                 |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                         |
| Species Reactivity         | Human                                                                                                                                                                                                                                                               |
| Product Description        | The GENLISA Human Oculomedin (OCLM) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use<br>Standardisation and High Reproducibility Lot to Lot Consistency<br>Accuracy and Precision Validated against seven points for a... |
| Range                      | 15.625 - 1000 pg/ml                                                                                                                                                                                                                                                 |
| Sensitivity                | 9.375 pg/ml                                                                                                                                                                                                                                                         |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                                  |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                          |