

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Nicotinamide Phosphoribosyltransferase (NAMPT) ELISA KBS-KBH1064

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Nicotinamide Phosphoribosyltransferase (NAMPT) ELISA                                                                                                                                                                                            |
| Biozol Catalog Number      | KBS-KBH1064                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | KBH1064                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | KBS-KBH1064-1                                                                                                                                                                                                                                                 |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Human                                                                                                                                                                                                                                                         |
| Product Description        | The GENLISA Human Nicotinamide Phosphoribosyltransferase (NAMPT) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Precision Validat... |
| Range                      | 0.315 - 20 ng/ml                                                                                                                                                                                                                                              |
| Sensitivity                | 0.15 ng/ml                                                                                                                                                                                                                                                    |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                            |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |