

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

AVB Sepharose High Performance CRM-AV-341C

Article Name	AVB Sepharose High Performance
Biozol Catalog Number	CRM-AV-341C
Supplier Catalog Number	AV-341C
Alternative Catalog Number	CRM-AV-341C-25
Manufacturer	Creative BioMart
Category	Sonstiges
Application Notes	<p>When using AVB Sepharose High Performance the AAV can be applied directly from clarified AAV vector cell lysate. Conventional buffers (e.g., PBS, Tris, citrate) may be used for loading, washing, and elution. Virus binds to the column at around neutral pH and is typically eluted by lowering the pH, for example in the range of pH 2 to 5. Since AAV is sensitive to highly acidic conditions, it is important to minimize the exposure to low pH during elution. Therefore, collected elution fractions should be neutralized immediately.</p>