

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Canine APOE (Apolipoprotein E) ELISA Kit EKL-ECA27286

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Canine APOE (Apolipoprotein E) ELISA Kit                                                                                                                                                                                                                      |
| Biozol Catalog Number      | EKL-ECA27286                                                                                                                                                                                                                                                  |
| Supplier Catalog Number    | ECa27286                                                                                                                                                                                                                                                      |
| Alternative Catalog Number | EKL-ECA27286-96T                                                                                                                                                                                                                                              |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Canine                                                                                                                                                                                                                                                        |
| Product Description        | This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Canine APOE. Standards or samples are added to the micro ELISA plate wells and combined with the specific an... |
| Range                      | 1.56-100ng/mL                                                                                                                                                                                                                                                 |
| Sensitivity                | 0.65ng/mL                                                                                                                                                                                                                                                     |
| Samples                    | Serum, plasma and other biological fluids                                                                                                                                                                                                                     |
| Target                     | APOE                                                                                                                                                                                                                                                          |