

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Porcine AP (Aprotinin) ELISA Kit EKL-EP22539

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Porcine AP (Aprotinin) ELISA Kit                                                                                                                                                                                                                              |
| Biozol Catalog Number      | EKL-EP22539                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | EP22539                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | EKL-EP22539-96T                                                                                                                                                                                                                                               |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Porcine                                                                                                                                                                                                                                                       |
| Product Description        | This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Porcine AP. During the reaction, Porcine AP in the sample or standard competes with a fixed amount of Porcine AP on the solid pha... |
| Range                      | 1.23-100ng/mL                                                                                                                                                                                                                                                 |
| Sensitivity                | 0.53ng/mL                                                                                                                                                                                                                                                     |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | AP                                                                                                                                                                                                                                                            |