

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-Human IgG-coupled Magnetic Beads (recommended for MPCLIA), Multiple ABS-MPC-A004-100MG

Article Name	Anti-Human IgG-coupled Magnetic Beads (recommended for MPCLIA), Multiple
Biozol Catalog Number	ABS-MPC-A004-100MG
Supplier Catalog Number	MPC-A004-100mg
Alternative Catalog Number	ABS-MPC-A004-100MG
Manufacturer	AcroBiosystems
Host	Multiple
Category	Sonstiges
Species Reactivity	Mouse
Product Description	AntiHuman IgGcoupled Magnetic Beads (recommended for MPCLIA) is expressed from Mouse Hybridoma. - Sonstiges...
Form	Powder
Target	AHG
Application Notes	The Anti-Human IgG-coupled Magnetic Beads is used to capture the human IgG, it can combination with Acridinium ester markers in chemiluminescence technology, The Acridinium ester markers such as Streptavidin-Acridinium ester can capture the biotinylated proteins or molecules, this allows detection of antigen and antibody binding or antibody screening.