

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Protein L-coupled Magnetic Beads (recommended for MPCLIA), E. coli ABS-MPC-A008-100MG

Article Name	Protein L-coupled Magnetic Beads (recommended for MPCLIA), E. coli
Biozol Catalog Number	ABS-MPC-A008-100MG
Supplier Catalog Number	MPC-A008-100mg
Alternative Catalog Number	ABS-MPC-A008-100MG
Manufacturer	AcroBiosystems
Host	E. coli
Category	Sonstiges
Species Reactivity	Bacteria
Product Description	Protein Lcoupled Magnetic Beads (recommended for MPCLIA) is expressed from Peptostreptococcus magnus E. coli. - Sonstiges...
Form	Powder
Target	Protein L
Application Notes	The Protein L-coupled Magnetic Beads is used to capture the most antibodies, it can combination with Acridinium ester markers in chemiluminescence technology. The Acridinium ester markers such as Streptavidin-Acridinium ester can capture the biotinylated proteins or molecules, this allows detection of Fc receptors and antibody binding.