

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-Human IgM-coupled Magnetic Beads (recommended for MPCLIA) ABS-MPC-A010-10MG

Article Name	Anti-Human IgM-coupled Magnetic Beads (recommended for MPCLIA)
Biozol Catalog Number	ABS-MPC-A010-10MG
Supplier Catalog Number	MPC-A010-10mg
Alternative Catalog Number	ABS-MPC-A010-10MG
Manufacturer	AcroBiosystems
Host	Human
Category	Sonstiges
Species Reactivity	Human
Product Description	AntiHuman IgMcoupled Magnetic Beads (recommended for MPCLIA) is expressed from Human HEK293. - Sonstiges...
Form	Powder
Target	Anti-Human IgM
Application Notes	The Anti-Human IgM-coupled Magnetic Beads is used to capture the human IgM, it can combination with Acridinium ester markers in chemiluminescence technology, The Acridinium ester markers such as Streptavidin-Acridinium ester can capture the biotinylated proteins or molecules, this allows detection of antigen and antibody binding or antibody screening.