

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### 10\* Capping Reaction Buffer ABN-U0855

|                            |                                                                                                                                    |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | 10* Capping Reaction Buffer                                                                                                        |
| Biozol Catalog Number      | ABN-U0855                                                                                                                          |
| Supplier Catalog Number    | U0855                                                                                                                              |
| Alternative Catalog Number | ABN-U0855-1                                                                                                                        |
| Manufacturer               | Abnova                                                                                                                             |
| Category                   | Sonstiges                                                                                                                          |
| Product Description        | 10*Capping Reaction Buffer is the companion reaction buffer for the Vaccinia Capping Enzyme and mRNA Cap 2-O-Methyltransferase.... |
| Buffer                     | 1*Capping Reaction Buffer is 50 mM Tris-HCl (pH 8.0 at 25C), 5 mM KCl, 1 mM DTT, 1 mM MgCl <sub>2</sub> .                          |
| Form                       | Liquid                                                                                                                             |
| Application Dilute         | The optimal working dilution should be determined by the end user.                                                                 |