

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### ELISA Kit for Homovanillic Acid (HVA) USC-CED717GE

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| Article Name               | ELISA Kit for Homovanillic Acid (HVA)                                                                |
| Biozol Catalog Number      | USC-CED717GE                                                                                         |
| Supplier Catalog Number    | CED717Ge                                                                                             |
| Alternative Catalog Number | USC-CED717GE-48, USC-CED717GE-96                                                                     |
| Manufacturer               | Cloud-Clone                                                                                          |
| Category                   | Kits/Assays                                                                                          |
| Species Reactivity         | All                                                                                                  |
| Product Description        | 48T, 96T, 96T*5, 96T*10, 96T*100...                                                                  |
| Range                      | 2.47-200ng/mL                                                                                        |
| Sensitivity                | The minimum detectable dose of this kit is typically less than 0.95ng/mL                             |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Application Notes          | Enzyme-linked immunosorbent assay for Antigen Detection.                                             |