

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Melanoma Marker Antibody Cocktail (MART-1 & Tyrosinase & gp100), Clone: [DT101+ BC199 + T311 + HMB45], Mouse, Monoclonal NSJ-V8328-20UG**

Article Name	Melanoma Marker Antibody Cocktail (MART-1 & Tyrosinase & gp100), Clone: [DT101+ BC199 + T311 + HMB45], Mouse, Monoclonal
Biozol Catalog Number	NSJ-V8328-20UG
Supplier Catalog Number	V8328-20UG
Alternative Catalog Number	NSJ-V8328-20UG
Manufacturer	NSJ Bioreagents
Host	Mouse
Category	Antikörper
Application	IHC-P
Species Reactivity	Human
Immunogen	Recombinant hMART-1 protein (DT101 & BC199), extract of pigmented melanoma metastases from lymph nodes (HMB45) and recombinant tyrosinase protein (T311) were used as the immunogen for this Melanoma Marker antibody cocktail.
Product Description	This antibody cocktail recognizes three melanoma-specific proteins, which include MART-1, Tyrosinase and gp100. MART-1 is a newly identified melanocyte differentiation antigen recognized by autologous cytotoxic T lymphocytes. Tyrosinase is one of the...
Clonality	Monoclonal
Clone Designation	[DT101+ BC199 + T311 + HMB45]
Purity	Protein G affinity chromatography

Form	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Antibody Type	Primary Antibody
Application Dilute	Immunohistochemistry (FFPE): 1-2ug/ml
Application Notes	Optimal dilution of the Melanoma Marker antibody cocktail should be determined by the researcher.